







Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931

Accessories of Slit lamp

Product name	Model#	Photo	Remark
Contact lens	YZ15		Can be used with any brand slit lamp.
Gonioscope- single lens	YZ14		Can be used with any brand slit lamp.
Tetrahedron Lens			Can be used with any brand slit lamp.
Trihedron lens	YZ13		Can be used with any brand slit lamp.
Aspheric Lenses 90D			Can be used with any brand slit lamp.
Aspheric Lenses 78D			Can be used with any brand slit lamp.
Aspheric Lenses 20D			Can be used with any brand slit lamp.

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931

Accessories of Slit lamp

Product name	Model#	Photo	Remark
Cable(connecting button and DC)			For after-service use. Support for S350DC and S350 Image System slit lamp.
LED module(for zeiss type)			3v/3w Support for S280 series slit lamp.
LED module(for tower type)			3v/3w Support for S350 series and S390H/L slit lamp.
Diffusion mirror(for zeiss type)			For after-service use. Support for S280 series slit lamp.
Diffusion mirror(for tower type)			For after-service use. Support for S350 series and S390H/L slit lamp.
Dustproof cover			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Movement plate			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Rail cover			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931









Accessories of Slit lamp

Product name	Model#	Photo	Remark
Rail			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Gear wheel			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Pin			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Lamp shade			For after-service use. Support for S350 series and S390H/L slit lamp.
Forehead band			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Filter			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Microscope complete set(5 magnifiers)			It includes eyepiece component, small and big objective lenses and 5-magnifier drum.
Microscope complete set(3 magnifiers)			It includes eyepiece component, small and big objective lenses and 3-magnifier drum.
Power cable			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp. 3 of 7

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931








Accessories of Slit lamp

Product name	Model#	Photo	Remark
Power supply			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Shield glass			For after-service use. Support for S350 series and S390H/L slit lamp.
Switch			Support for the electric table(AT-16) of standard slit lamp.
Fuse			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Chin paper			400pcs/set
Mirror lens			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
Eyepiece			Support for S350 series, S280 series and S390H/L slit lamp.
ORSAM lamp	HLX6425 1		6V/20W, support for S350 series and S280 series slit lamp. 4 of 7

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931








Accessories of Slit lamp

Product name	Model#	Photo	Remark
Black metal connecting board			It's matched under optometry table.
Aluminium plate			It's matched on optometry table.
Goldmann T type tonometer adaptor			supported by slit lamp models: S350 S350s
Goldmann R type tonometer adaptor			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.
870 tonometer adaptor			supported by slit lamp models: S280
Baclite background illumination model			Improve the image quality of digital images, especially in the narrow crack of light inspection. (supported by slit lamp model: S280 S350 S350s)
Teaching tube			other doctors, students, nurses or technical personnel can be involved in the inspection process through teaching tube. (supported by slit lamp model: S280 S350 S350s)

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931





Accessories of Slit lamp

Product name	Model#	Photo	Remark
Yellow filter module			Improve contrast of image when fluorescence staining, can be permanently installed in the microscope,use it or not by sliding switch.(supported by slit lamp model: S280 S350 S350s)
Polarizer Module			Used to reduce the corneal surface reflected light,Can be permanently installed in the slit lamp microscope, use it or not by sliding switch.(supported by slit lamp model: S280 S350 S350s.)
Retical(Measuring eyepiece)			Can be installed in the eyepiece, upgrade the common eyepiece to measurement eyepiece. (supported by slit lamp model: S280 S350 S350s)
Elbow Supporting Module			Help doctors alleviate elbow fatigue,4pcs/set, through the addition and subtraction number of pieces by different working height.
Foot-controlled shooting model			Support for S350 Image System and S390H/L slit lamp.
Amber filter			By improving the color depth to obtain a more clear and higher contrast images, can be permanently installed in the lighting components of slit lampit needs to be installed and commissioned by factory, please consult sales staff before ordering. (supported by slit lamp models: S350 S350s)
Yellow filter			For after-service use. Support for S350 series, S280 series and S390H/L slit lamp.

Shanghai Jiahe Optical Co., Ltd

Website:www.jhoptical.com Tel:0086-21-60481633 Fax:0086-21-33321931

Accessories of Slit lamp

Product name	Model#	Photo	Remark
Image process software			Support for S350DC, S350 Image System and S390H/L slit lamp.
CCD adaptor			C-type interface, including focus tuning and aperture adjustment function. (supported by slit lamp model: S280 S350 S350s)
Digital camera adaptor			Canon EF mount. (supported by slit lamp model: S280 S350 S350s)
Beam splitter			As shown, there are two types of splitter structures to choose from, please confirm your slit lamp microscope models before ordering. (supported by slit lamp model: S280 S350 S350s)